


BM02



30 Sh., 6 mm

25 Sh., 15 mm

For children preventing insoles

BM03



35 Sh., 15 mm

25 Sh., 15 mm

For elastic sport insoles

BM04



25 Sh., 15 mm

25 Sh., 15 mm

For patients with diabetes and rheumatism

BM06



60 Sh., 6 mm

25 Sh., 15 mm

For daily orthopedic adult insoles

BM07




60 Sh., 6 mm

35 Sh., 15 mm

For daily orthopedic adult insoles

BM08



50 Sh., 2 mm

15 Sh., 6,4 mm

35 Sh., 15 mm

For elastic sport insoles with extra shock absorber layer

BM09



45 Sh., 40 mm

For patients with diabetes and rheumatism and cerebral paralysis

BM10



35 Sh., 40 mm

For patients with diabetes and rheumatism and cerebral paralysis

BM11



55 Sh., 40 mm

For daily orthopedic adult men insoles with weight over 75 kg

BM13



35/55 Sh., 40 mm

For daily orthopedic adult with extra midfoot suport

BM14



30 Sh., 2 mm

25 Sh., 15 mm

For children preventing insoles

BM15



30 Sh., 2 mm

55 Sh., 40 mm

For daily orthopedic adult men insoles with weight over 75 kg with schock absorbing layer