









ORTOPEDIJOS TECHNIKA

MATERIALS

2023-1





Leathers	1
Textiles	8
Linings	10
Paddings	12
Stiffenings	13
Fastenings	15
Shoelaces	17
Leather welts	18
Welts	19
Midsole materials	20
CNC midsoles	21
Boardings	22
Soles	23
Molded soles	29
Wooden heels	30
Vacuum insoles	31
Milled insoles	32
Insoles top covers	33



LILY

Thickness
1.0-1.2 mm



SETA

Thickness
1.3-1.5 mm



TULIP

Thickness
1.1-1.3 mm



ART

Thickness
1.1-1.2 mm



MOUSSE


Thickness
1.0-1.1 mm





SNAKY

Thickness
0.9-1.1 mm




VELVET

Thickness
0.9-1.0 mm



VOLCANO

Thickness
0.9-1.1 mm



PURESOFT

Thickness
1.1-1.3 mm






LUMIÈRE - LIPSTICK

Thickness
0.9-1.1 mm



CRONOS

Thickness
0.9-1.1 mm




VANESSA

Thickness
0.8-1.0 mm




COCCO OYSTER

Thickness
0.9-1.0 mm




SOUPLESSE

Thickness
0.9-1.0 mm



BOMBAY


Thickness
0.9-1.0 mm





JAIPUR PEARL

Thickness
0.9-1.1 mm




VEDETTE

Thickness
1.1-1.3 mm




KIDEPO

Thickness
0.8-0.9 mm



MEDICI

Thickness
1.2-1.4 mm



BOXLIGHT

Thickness
1.0-1.2 mm




COLORADO

Thickness
1.1-1.3 mm




IMPRUNETA

Thickness
1.4-1.6 mm



SPORTLUX

Thickness
1.3-1.5 mm



HUNTER

Thickness
1.6-1.8 mm



**SURVIVAL
HYDRO**

Thickness
1.4-1.6 mm



ESPRIT

Thickness
1.2-1.4 mm



TIFFANY
HYDRO

Thickness
1.3-1.5 mm



ONTARIO
HYDRO

Thickness
1.5-1.7 mm




PAMPA
HYDRO

Thickness
1.7-1.9 mm



HEAVY
DUTY

Thickness
2.0-2.2 mm





SPECIAL MATERIALS






TX21
White



TX20
Light grey




TX22
Gris grey



TX17
Green with
padding




TX10
Red




TX11
Azur




TX23
Azul blue




TX24
Marine




TX07
Black




TX27
Chaki with
padding




TX25
Jeans with
padding




TX26
Black with
padding




TX18
Elastic



TX19
Tricot negro



TX28
Recycled
black



TX29
Stretchable
black



TX30
Camouflage



TX31
Cordura black




Recycled material



TX32
Cordura brown




TX33
Cordura blue




TX34
Cordura chaki



TX35
Cordura grey



TX37
Cordura granate



LN02

Calflining cream,
0.7-0.9 mm



LN45

Natura lining
burro,
1.0-1.1 mm




LN26

Lamnappa
naturel,
0.7-0.9 mm



LN22

Calflining camel,
0.9-1.1 mm



LN24

Calflining nero,
0.7-0.9 mm



LN46

Recycled
camel,
1 mm



LN47

DuraFit PLUS,
1.1 mm




LN43

Foam laminated
PE smooth
suede-like finish
beige,
4 mm




LN38

Foam laminated
PE smooth
suede-like finish
chocolate,
4 mm




LN39

Foam laminated
PE smooth
suede-like finish
nero,
4 mm




LN06

Foam laminated
PE flisse black,
4 mm



LN44

Foam laminated
PE mali velour
dark grey,
4 mm



LN48

Foam laminated
black



Recycled material




Chrom free




LN42
Puratex,
3.2 mm




LN28
Natural fur
off white



LN30
Natural fur
brown



LN27
Natural fur
black



LN36
Wool grey lining



Waterproof

LATEX

SFM 01-09
8 mm

SFM 02-09
3 mm

SFM 03-09
5 mm



PU FOAM

SFM 04-81
3 mm

SFM 05-81
5 mm

SFM 06-81
8 mm




ZEG

SFM 07-81
3 mm


SFM 08-81
5 mm



TOE CAP



Interflex



MAX



Erkoflex



TEPP2




Polypropylene transparent



Simona




Pre-preg




Metal toe cap




HEEL COUNTER



Interflex




Erkoflex



Leather



MAX



Rhenoflex



TEPP2



see page 14 for more heel counter materials

Rubberflex




Polypropylene transparent




Simona



Pre-preg




BRANDSOLE



Ibiflex




Thermoflex




Leather



Pre-preg



Carbon



Rubberflex




Eva
35 Sh



EYELET




HOOK






TEXTILE HOOK




STUD



BUCKLE



QUICK LASING SYSTEM



CORD LOCKS



**COTTON
ROUND**
Length 75 cm



White

Grey

Yellow



Blue

Brown

Black

**POLYESTER
ROUND**
Length
60, 75, 90, 120,
150, 180, 200 cm



Sand

**COTTON
MULTICOLOR**
Length
90, 120 cm



Blue with white

Grey with red

**COTTON
FLAT**
Length
75, 90 cm



Sand

Red

Blue



Brown

White (+ extra 120 cm)

Black (+ extra 120cm)

**COTTON
WAXED**
Length 75 cm



Black

**ELASTIC
ROUND**
Length 60 cm



Black

White with black



7410021

7410020



7411210

7411214



7411241

7411240



7411290

EVA

RUBBER



1-09



1-82



1-81



3-19



3-17



3-41



3-46



3-56



3-78



3-81



4-23



4-25



4-26



6-09



6-17



6-46



6-81

RUBBER



24-46




24-81



26-19



26-17



26-46



26-81



31-17



31-46



31-81



32-46



32-81



36-17



36-19



36-81



39-46



39-81



42-81



44-46




44-81




EVA
09 White



EVA
19 Light beige



EVA
17 Dark beige




EVA
41 Taupe



EVA
56 Grey



EVA
46 Brown



EVA
81 Black



Benefit of CNC soles:

- no visible gluing seams
- lightness

Usage of CNC soles:

- for shoes with roller
- sport type and casual shoes

Available colours:

- Black 81

The title of material in Odoo system:

CNC A4H (Flat + Boarding)



CNC A4R (Flat + Welt)



CNC A6R (Formed + Welt)

Available hardness:

- EVA 45 Sh
- EVA 55 Sh

How to order in Odoo?

1. Select Lasting type „Cemented without welt“.
2. Select outer sole type from drop down list.
3. Specify requested heel height.
4. Specify requested hardness.
5. Standard roller thickness 5 mm, if other thickness requested please specify.
6. In comment field write a type of requested CNC sole.


RUBBER
1,8 mm



BR31-87

BR30-81

EVA
2 mm



BR08-19

BR07-17

BR06-41

BR05-46

BR04-81

TPU
1,2 mm



BR23-09/89

BR24-09/81

BR25-56/81

DIAMOND
1,8 mm



BR22-09

BR28-81

BR29-09

LAVERO
3 mm



BR20-09

BR18-19

BR21-05



BR19-17

BR15-56

BR17-41

BR16-46

BR14-81

STONE
Recycled
3 mm



BR27-81

SPORT




BR26-81




Recycled material


Super New Flex
35 Sh



Blue Clasik 701
50 Sh





Vibram 8870
45 Sh




Astro-Star
50 Sh



Swallow
45 Sh




Vibram Mandorlo
45 Sh



Astrolight
50 Sh




Zephir VIP Classik
50 Sh



Pluriball
New Flex
50 Sh




Astrol
65 Sh



Stone
Recycled
50 Sh



Autoband
65 Sh




Topy
Cellolux
55 Sh




Hexagon
65 Sh




Crepina
Soft
60 Sh




Contol
65 Sh



Vibram
Claw
60 Sh



Nylon Topy



Resistant to oil and acid



Recycled material




Vibram 446



Vibram 1220




Vibram 1319



Vibram 2002



Vibram 2609



Vibram 2644



Vibram 2667



Vibram 2668



Vibram 2690





Vibram 528



Vibram 2603



Vibram 2655



Vibram 2678



Vibram 2810



Flensburg



Flex



Grip Line 505



Kiel



Resistant to oil and acid



Antistatic



Lubeck



81

Mammoth



87

46

81

Sportflex



09

46

81

5001



87

46

81

Venus



87

81

Ventosa



09

87

78

Carrion



81

Vibram 586



09

81

Vibram 9105



09

81



Resistant to oil and acid



Antistatic

Roccia Block



Sole leather



Anti-slip Sole - 52



Anti-slip Heel - 52



Anti-slip Half sole



Shell PU1662

Size
40-45




Nadia

Size
36-41



Nel

Size
37-41





Nr. 2



Nr. 6



LEFA



← Heel

← Heel tap


Cork

3 mm
5 mm
10 mm
15 mm



Microcork

60 Sh
5 mm
10 mm



Nora

30 Sh
4 mm
6 mm




Multiform

30 Sh
1.6 mm



Multiform

35 Sh
perforated
3.5 mm




Multiform

35 Sh
perforated
3.5 mm



Multiform

25 Sh
15 mm



Multiform

35 Sh
5 mm
7 mm
10 mm
15 mm



Multiform

16 Sh
perforated
6 mm



PPT

1.6 mm



PPT

3.2 mm




PPT

6.4 mm



Poron

1.6 mm



Poron

3.2 mm



Plastazote

3 mm
8 mm
12 mm



New flex

6 mm



Multiform

35 Sh
10 mm



Microcork

35 Sh
5 mm
10 mm




BM02



For children preventing insoles

BM03



For elastic sport insoles

BM04



For patients with diabetes and rheumatism

BM06




For daily orthopedic adult insoles

BM07



For daily orthopedic adult insoles

BM08




For elastic sport insoles with extra shock absorber layer

BM09



For patients with diabetes and rheumatism and cerebral paralysis

BM10




For patients with diabetes and rheumatism and cerebral paralysis

BM11



For daily orthopedic adult men insoles with weight over 75 kg

BM13



For daily orthopedic adult with extra midfoot suport

BM14



For children preventing insoles

BM15



For daily orthopedic adult men insoles with weight over 75 kg with schock absorbing layer

LN32
Horse leather,
0.9-1.1 mm



LN02
Calflining cream,
0.7-0.9 mm



LN45
Natura lining burro,
1.0-1.1 mm



LN26
Lamnappa naturel,
0.7-0.9 mm



LN22
Calflining camel,
0.9-1.1 mm



LN24
Calflining nero,
0.7-0.9 mm



LN18
Suetex beige,
0.7 mm



LN17
Suetex grey,
0.7 mm




LN19
Micro suetex
black,
0.6 mm



LN33
ON STEAM
perforated
beige,
0.8 mm



LN34
ON STEAM
perforated
black,
0.8 mm




LN46
Recycled
camel,
1 mm



LN47
DuraFit PLUS,
1.1 mm



LN49
Medi-fiber blue,
0.5 mm






LN28

Natural fur
off white



LN30

Natural fur
brown



LN27

Natural fur
black



LN44

Foam laminated
PE mali velour
dark grey,
4 mm





ORTOPEDIJOS TECHNIKA

Krėvės pr. 129, LT-50312 Kaunas, Lithuania
E-mail: info@orthopedic-pro.com

Phone +370 37 312 686
Mobile +370 665 28 542

www.orthopedic-pro.com

