



YOU ASK



ORTOPEDIJOS TECHNIKA



WE MAKE



THEY  
WALK



*Happy customers worldwide!*

# MANUFACTURING TRADITIONS SINCE 1945!

- Leading market position.
- More than **600 experienced employees**.
- Ortopedijos Technika exports their products to **24 countries**.
- ISO 9001** quality management system since 2001.
- Compliance with requirements of **Regulation (EU) 2017/745 on medical devices (MDR)** and **FDA**.



## *Your challenges - our solutions!*

We offer not only orthopedic products, production services, but first of all - individual decisions for each client.

### **The creation and development of our products is based on:**

- teamwork of our therapists, orthopedic technicians and production technologists;
- feedback from our patients;
- ideas or requests from our customers.



ORTOPEDIJOS TECHNIKA



Share your challenges with us and we'll offer you the best solution. You can do this by calling or writing to us and in person at the international trade fairs and congresses *OT World Leipzig* and *Orthopädie Schuh Technik Cologne*.



ORTOPEDIJOS TECHNIKA

# Pre-preg technologies

Maximum strength and durability;

Strong and lightweight solution to provide constant support to the foot;

Sole can be easily adapted to different shoes.


## Custom made products:

ankle-foot orthosis (AFO);

knee orthosis (KO);

knee-ankle-foot orthosis (KAFO);

with or without joints.



“WATER PRINT technology provides an opportunity to personalize the product and to create BEAUTY and STYLE”



## Standard Pre-preg orthosis

**4 models of standard AFO (anterior, posterior lateral and medial)**

- CE mark;
- made only from carbon fibre or from carbon and glass fibre;
- with or without soft parts;
- with faster or slower rocker sole;
- convenient packaging solution.




Orthosis could be with or without trimmable sole part.



Pre-preg footplate could be manufactured serial as well as custom made.

Development of customers' models.  
Private label.







ORTOPEDIJOS TECHNIKA

# Orthopedic footwear

Our company has installed modern CAD-CAM technologies involving design and milling of shoe lasts and insoles, construction of uppers and cutting out of footwear elements.


## OPERATIONS OF TECHNOLOGICAL PRODUCTION PROCESS:

- **Ordering:** via on-line ordering system. Patient data could be provided: by pedilen imprint / blueprint, plaster negatives, 3D scanned foot data;
- **Design and milling of shoe lasts;**
- **Modelling and milling of insoles;**
- **Production of test shoes;**
- **Pattern engineering;**
- **Cutting and sewing;**
- **Possibility to order the uppers without lasts sending;**
- **Lasting** (shaping the upper around the last);
- **Bottoming and finishing.**

Set of all materials – available by request.

**The quality** of our services is guaranteed by our quality management system.

Old manufacturing traditions and manual work are combined with extraordinary knowledge and modern technologies in order to offer the best result for You.



*Made for comfort, stability, support and still elegant*



ORTOPEDIJOS TECHNIKA

- professional solutions for complex cases;
- loose lining shoes;
- different types of heel and toe stiffenings;
- integrated orthoses (from carbon fiber and other materials);
- temporary post – operative shoes;
- orthopedic safety shoes;
- orthopedic shoes for sports.



# Footwear insoles

Various kinds of serial or custom made, milled or vacuum formed footwear insoles:

- for people with diabetes and rheumatism;**
- for children;**
- for athletes;**
- for each one of us.**


More than 15 standard models.

Development according the idea of the customer or sample.

Private label.

Set of different insole materials – available upon request.

**FIXALTO®**  
**StepFIX™**











ORTOPEDIJOS TECHNIKA

We offer You precise and quick design and production of insoles by using our modern CAD-CAM equipment.

**The following data may be used for the production of insoles:**

- digital data foot scan (via-e-mail);
- digital data of pedilen imprint, blueprint, plaster cast or shoe last (via e-mail);
- pedilen imprint, blueprint, plaster cast or shoe last (via post).

**Please Note:** Special customer's requirements for the insoles finishing could be implemented.



## Braces

### Knitted or/and sewed from elastic materials

Rehabilitation and sport braces which are used in cases of various diseases, traumas or after surgical operations are necessary for stabilizing different parts of the body in the right position and decreasing physical load on them.

Our braces provide the following characteristics: swelling control, proper compression, support and decreasing pain, staying mobile and stable at the same time.

**Note.** If you can't choose products from our catalogues we are eager to start manufacturing of new products. Production under private label is possible as well.

**RehaFIX®****SpecialFIX™****SoftFIX®****BasicFIX™**



### Leading market position

*Ortopedijos Technika* is an absolute market leader in the Lithuanian orthopedic market holding more than 37% market share in the country.




### Full range of products and services

The product portfolio of *Ortopedijos Technika* includes all sorts of orthopedic articles (prosthetic appliances, footwear, knitted products, etc.). *Ortopedijos Technika* also provides various orthopedic services to their clients.



### Experienced employees

Being the pioneer in the market *Ortopedijos Technika* became a key expert in the field of orthopedics in Lithuania. Employees of *Ortopedijos Technika* have an extensive long term experience.



### Competitive pricing

Economies of scale, high quality products and efficient production organization enable *Ortopedijos Technika* to offer its products and services under very competitive terms.



### Potential in foreign markets

*Ortopedijos Technika* has established business relationship with many foreign companies: 79 partners from 24 different countries.



ORTOPEDIJOS **TECHNIKA**

Krėvės pr. 129,  
LT-50312 Kaunas, Lithuania

E-mail: [info@orthopedic-pro.com](mailto:info@orthopedic-pro.com)  
[www.orthopedic-pro.com](http://www.orthopedic-pro.com)

Tel. +370 37 312 686  
Mobile +370 665 28 542



*Let's CREATE TOGETHER!*

